



2017

**Stealth Products**  
*Solutions For Life*®



# AL BR 10 L A R

Shell type

Back rest

Chair width

Height

Contour

Mounting

**CF - CF Series**

**AL - AL Series**

**AFT - AFT Series**

**GU - G Series**

**Note:**

Shell Assembly Weight

AL Series - 1.5lb / 680gm

AFT Series - 3.0lb / 1390gm

CF Series - .820lb / 370gm

- above weights are of a 16" STD

- all weights approx.

**10" / 25cm**

**11" / 28cm**

**12" / 30.5cm**

**13" / 33cm**

**14" / 35.5cm**

**15" / 38cm**

**16" / 40.5cm**

**17" / 43cm**

**L**  
low  
(10"/25cm)

**M**  
medium  
(13"/33cm)

**S**  
standard  
(16"/40.5cm)

**T**  
tall  
(18"/45.7cm)

**X**  
x-tall  
(20"/51cm)

**A**  
active  
(2.5"/63.5mm)

**D**  
deep  
(4"/100mm)

**6**  
extra deep  
(6"/150mm)

**F**  
Fixed Elite  
(.60lb / 272gm)

**2**  
2-Point Pro  
(1.25lb/566gm)

**4**  
4-Point Pro  
(1.25lb/566gm)

**R**  
Quick  
Release  
(1.30lb / 589gm)

**W**  
Width  
Adjustable Elite  
(.63lb / 285gm)

**B**  
Basic Release  
(1.50lb / 680gm)



# CF Series Back Supports

The CF Series Back Support utilizes a high quality carbon fiber shell that is not only ultra lightweight, but also highly durable and attractive.

#### Fits Chair Widths

12" (30.5cm) - 19" (48cm)

#### Support Heights

Low: 10" (25cm)

Medium: 13" (33cm)

Standard: 16" (41cm)

#### Contour Depths

Active: 2.5" (63.5mm)

Deep 4" (100mm)

#### Weight

Weight: .820lb (370 Grams)

(approx. weight of a 16" / 40.5cm standard height back)

#### NOTE:

Custom heights available





User Inspired Designs

# CF Series

## Low Height (10"/25cm), Positioning Back Supports

### Active Contour Depth

#### Performance Mounting Hardware Options

Width of Wheelchair	Fixed ELITE	Width Adjustable ELITE	Quick Release	2-Point PRO	4-Point PRO	MSRP
12" / 30.5cm	CFBR12LAF	CFBR12LAW	CFBR12LAR	CFBR12LA2	CFBR12LA4	\$745.00
13" / 33cm	CFBR13LAF	CFBR13LAW	CFBR13LAR	CFBR13LA2	CFBR13LA4	
14" / 35.5cm	CFBR14LAF	CFBR14LAW	CFBR14LAR	CFBR14LA2	CFBR14LA4	
15" / 38cm	CFBR15LAF	CFBR15LAW	CFBR15LAR	CFBR15LA2	CFBR15LA4	
16" / 40.6cm	CFBR16LAF	CFBR16LAW	CFBR16LAR	CFBR16LA2	CFBR16LA4	
17" / 43cm	CFBR17LAF	CFBR17LAW	CFBR17LAR	CFBR17LA2	CFBR17LA4	
18" / 45.7cm	CFBR18LAF	CFBR18LAW	CFBR18LAR	CFBR18LA2	CFBR18LA4	
19" / 48.2cm	n/a	CFBR19LAW	CFBR19LAR	CFBR19LA2	CFBR19LA4	

### Deep Contour Depth, (4"/100mm), Coded K0108

Width of Wheelchair	Fixed ELITE	Width Adjustable ELITE	Quick Release	2-Point PRO	4-Point PRO	MSRP
12" / 30.5cm	CFBR12LDF	CFBR12LDW	CFBR12LDR	CFBR12LD2	CFBR12LD4	\$745.00
13" / 33cm	CFBR13LDF	CFBR13LDW	CFBR13LDR	CFBR13LD2	CFBR13LD4	
14" / 35.5cm	CFBR14LDF	CFBR14LDW	CFBR14LDR	CFBR14LD2	CFBR14LD4	
15" / 38cm	CFBR15LDF	CFBR15LDW	CFBR15LDR	CFBR15LD2	CFBR15LD4	
16" / 40.6cm	CFBR16LDF	CFBR16LDW	CFBR16LDR	CFBR16LD2	CFBR16LD4	
17" / 43cm	CFBR17LDF	CFBR17LDW	CFBR17LDR	CFBR17LD2	CFBR17LD4	
18" / 45.7cm	CFBR18LDF	CFBR18LDW	CFBR18LDR	CFBR18LD2	CFBR18LD4	
19" / 48.2cm	n/a	CFBR19LDW	CFBR19LDR	CFBR19LD2	CFBR19LD4	



**NOTE:**

- 1) Each back support includes: support shell, mounting hardware, contoured foam cushion, cover, and operation manual.
- 2) 19" CF Series backs utilize the same width shell supports as 18" CF Series backs



# CF Series

## Medium Height(13" / 33cm), Positioning Back Support Active Contour Depth(2.5" / 63.5mm), Coded K0108

### Performance Mounting Hardware Options

Width of Wheelchair	Fixed ELITE	Width Adjustable ELITE	Quick Release	2-point PRO	4-Point PRO	MSRP
12" / 30.5cm	CFBR12MAF	CFBR12MAW	CFBR12MAR	CFBR12MA2	CFBR12MA4	\$815.00
13" / 33cm	CFBR13MAF	CFBR13MAW	CFBR13MAR	CFBR13MA2	CFBR13MA4	
14" / 35.5cm	CFBR14MAF	CFBR14MAW	CFBR14MAR	CFBR14MA2	CFBR14MA4	
15" / 38cm	CFBR15MAF	CFBR15MAW	CFBR15MAR	CFBR15MA2	CFBR15MA4	
16" / 40.6cm	CFBR16MAF	CFBR16MAW	CFBR16MAR	CFBR16MA2	CFBR16MA4	
17" / 43cm	CFBR17MAF	CFBR17MAW	CFBR17MAR	CFBR17MA2	CFBR17MA4	
18" / 45.7cm	CFBR18MAF	CFBR18MAW	CFBR18MAR	CFBR18MA2	CFBR18MA4	
19" / 48.2cm	n/a	CFBR19MAW	CFBR19MAR	CFBR19MA2	CFBR19MA4	

## Deep Contour Depth(4" / 100mm), Coded K0108

### Performance Mounting Hardware Options

Width of Wheelchair	Fixed ELITE	Width Adjustable ELITE	Quick Release	2-point PRO	4-Point PRO	MSRP
12" / 30.5cm	CFBR12MDF	CFBR12MDW	CFBR12MDR	CFBR12MD2	CFBR12MD4	\$890.00
13" / 33cm	CFBR13MDF	CFBR13MDW	CFBR13MDR	CFBR13MD2	CFBR13MD4	
14" / 35.5cm	CFBR14MDF	CFBR14MDW	CFBR14MDR	CFBR14MD2	CFBR14MD4	
15" / 38cm	CFBR15MDF	CFBR15MDW	CFBR15MDR	CFBR15MD2	CFBR15MD4	
16" / 40.6cm	CFBR16MDF	CFBR16MDW	CFBR16MDR	CFBR16MD2	CFBR16MD4	
17" / 43cm	CFBR17MDF	CFBR17MDW	CFBR17MDR	CFBR17MD2	CFBR17MD4	
18" / 45.7cm	CFBR18MDF	CFBR18MDW	CFBR18MDR	CFBR18MD2	CFBR18MD4	
19" / 48.2cm	n/a	CFBR19MDW	CFBR19MDR	CFBR19MD2	CFBR19MD4	



**NOTE:**

- 1) Each back support includes: support shell, mounting hardware, contoured foam cushion, cover, and operation manual.
- 2) 19" CF Series backs utilize the same width shell supports as 18" CF Series backs

# CF Series

## Standard Height(16" / 41cm), Positioning Back Support

### Active Contour Depth(2.5"/63.5mm), Coded K0108

#### Mounting Hardware Options

Width of Wheelchair	Fixed ELITE	Width Adjustable ELITE	Quick Release	2-point PRO	4-Point PRO	MSRP
12" / 30.5cm	CFBR12AF	CFBR12SAW	CFBR12SAR	CFBR12SA2	CFBR12SA4	\$880.00
13" / 33cm	CFBR13AF	CFBR13SAW	CFBR13SAR	CFBR13SA2	CFBR13SA4	
14" / 35.5cm	CFBR14AF	CFBR14SAW	CFBR14SAR	CFBR14SA2	CFBR14SA4	
15" / 38cm	CFBR15AF	CFBR15SAW	CFBR15SAR	CFBR15SA2	CFBR15SA4	
16" / 40.6cm	CFBR16AF	CFBR16SAW	CFBR16SAR	CFBR16SA2	CFBR16SA4	
17" / 43cm	CFBR17AF	CFBR17SAW	CFBR17SAR	CFBR17SA2	CFBR17SA4	
18" / 45.7cm	CFBR18AF	CFBR18SAW	CFBR18SAR	CFBR18SA2	CFBR18SA4	
19" / 48.2cm	n/a	CFBR19SAW	CFBR19SAR	CFBR19SA2	CFBR19SA4	

### Deep Contour Depth(4"/100mm), Coded K0108

#### Performance Mounting Hardware Options

Width of Wheelchair	Fixed ELITE	Width Adjustable ELITE	Quick Release	2-point PRO	4-Point PRO	MSRP
12" / 30.5cm	CFBR12SDF	CFBR12SDW	CFBR12SDR	CFBR12SD2	CFBR12SD4	\$955.00
13" / 33cm	CFBR13SDF	CFBR13SDW	CFBR13SDR	CFBR13SD2	CFBR13SD4	
14" / 35.5cm	CFBR14SDF	CFBR14SDW	CFBR14SDR	CFBR14SD2	CFBR14SD4	
15" / 38cm	CFBR15SDF	CFBR15SDW	CFBR15SDR	CFBR15SD2	CFBR15SD4	
16" / 40.6cm	CFBR16SDF	CFBR16SDW	CFBR16SDR	CFBR16SD2	CFBR16SD4	
17" / 43cm	CFBR17SDF	CFBR17SDW	CFBR17SDR	CFBR17SD2	CFBR17SD4	
18" / 45.7cm	CFBR18SDF	CFBR18SDW	CFBR18SDR	CFBR18SD2	CFBR18SD4	
19" / 48.2cm	n/a	CFBR19SDW	CFBR19SDR	CFBR19SD2	CFBR19SD4	



**NOTE:**

- 1) Each back support includes: support shell, mounting hardware, contoured foam cushion, cover, and operation manual.
- 2) 19" CF Series backs utilize the same width shell supports as 18" CF Series backs

# AL Series Back Supports

**The AL Series Back Support provides the user with a lightweight, supportive backrest. Offered in a broad range of widths, heights and contours. Ideal for manual and power chair users alike.**

#### **Fits Chair Widths**

10" (25.4cm) - 20" (51cm)

#### **Support Heights**

Low: 10" (25cm)

Medium: 13" (33cm)

Standard: 16" (41cm)

Tall: 18" (43cm)

Extra Tall: 20" (51cm)

#### **Contour Depths**

Active: 2.5" (63.5mm)

Deep: 4" (100mm)

Extra Deep: 6" (150mm)

#### **Weight**

Weight: 1.5lbs (680 grams)

(approx. weight of a 16" / 40.5cm standard height back)

#### **HCPCS Coding**

Active - E2613

Deep - E2615

Extra Deep - E2615

#### **NOTE:**

Custom options available





User Inspired Designs

# AFT Series

## Tall Height(18"/46cm), Positioning Back Supports

### Active Contour Depth(2.5"/63.5mm), Code E 2613\*\*

#### Performance Mounting Hardware Options

#### Basic Mounting

Width of Wheelchair	2-point PRO	4-Point PRO	Quick Release	Width Adjustable ELITE	MSRP	Basic Release	MSRP
16" / 40.6cm	AFTBR16TA2	AFTBR16TA4	AFTBR16TAR	AFTBR16TAW	\$655.00	AFTBR16TAB	\$655.00
17" / 43cm	AFTBR17TA2	AFTBR17TA4	AFTBR17TAR	AFTBR17TAW		AFTBR17TAB	
18" / 45.7cm	AFTBR18TA2	AFTBR18TA4	AFTBR18TAR	AFTBR18TAW		AFTBR18TAB	
19" / 48.2cm	AFTBR19TA2	AFTBR19TA4	AFTBR19TAR	AFTBR19TAW		AFTBR19TAB	
20" / 51cm	AFTBR20TA2	AFTBR20TA4	AFTBR20TAR	AFTBR20TAW		AFTBR20TAB	

### Deep Contour Depth(4"/100mm), Code E2615\*

#### Performance Mounting Hardware Options

#### Basic Mounting


Width of Wheelchair	2-point PRO	4-Point PRO	Quick Release	Width Adjustable ELITE	MSRP	Basic Release	MSRP
16" / 40.6cm	AFTBR16TD2	AFTBR16TD4	AFTBR16TDR	AFTBR16T6W	\$730.00	AFTBR16TDB	\$730.00
17" / 43cm	AFTBR17TD2	AFTBR17TD4	AFTBR17TDR	AFTBR17T6W		AFTBR17TDB	
18" / 45.7cm	AFTBR18TD2	AFTBR18TD4	AFTBR18TDR	AFTBR18T6W		AFTBR18TDB	
19" / 48.2cm	AFTBR19TD2	AFTBR19TD4	AFTBR19TDR	AFTBR19T6W		AFTBR19TDB	
20" / 51cm	AFTBR20TD2	AFTBR20TD4	AFTBR20TDR	AFTBR20T6W		AFTBR20TDB	

### Extra Deep Contour Depth(6"/150mm), Code E2615

#### Performance Mounting Hardware Options

#### Basic Mounting

Width of Wheelchair	2-point PRO	4-Point PRO	Quick Release	Width Adjustable ELITE	MSRP	Basic Release	MSRP
16" / 40.6cm	AFTBR16T62	AFTBR16T64	AFTBR16T6R	AFTBR16T6W	\$805.00	AFTBR16T6B	\$805.00
17" / 43cm	AFTBR17T62	AFTBR17T64	AFTBR17T6R	AFTBR17T6W		AFTBR17T6B	
18" / 45.7cm	AFTBR18T62	AFTBR18T64	AFTBR18T6R	AFTBR18T6W		AFTBR18T6B	
19" / 48.2cm	AFTBR19T62	AFTBR19T64	AFTBR19T6R	AFTBR19T6W		AFTBR19T6B	
20" / 51cm	AFTBR20T62	AFTBR20T64	AFTBR20T6R	AFTBR20T6W		AFTBR20T6B	



**NOTE:**

1) Each back support includes: support shell, mounting hardware, contoured foam cushion, cover, and operation manual.





User Inspired Designs

# AL Series

## Low Height(10"/25cm). Positioning Back Supports

### Active Contour Depth(2.5"/63.5mm), Coded E2613

#### Performance Mounting Hardware Options

Width of Wheelchair	2-point PRO	4-Point PRO	Quick Release	Width Adjustable ELITE	MSRP
10" / 25.4cm	ALBR10LA2	ALBR10LA4	ALBR10LAR	ALBR10LAW	\$610.00
11" / 28cm	ALBR11LA2	ALBR11LA4	ALBR11LAR	ALBR11LAW	
12" / 30.5cm	ALBR12LA2	ALBR12LA4	ALBR12LAR	ALBR12LAW	
13" / 33cm	ALBR13LA2	ALBR13LA4	ALBR13LAR	ALBR13LAW	
14" / 35.5cm	ALBR14LA2	ALBR14LA4	ALBR14LAR	ALBR14LAW	
15" / 38cm	ALBR15LA2	ALBR15LA4	ALBR15LAR	ALBR15LAW	
16" / 40.6cm	ALBR16LA2	ALBR16LA4	ALBR16LAR	ALBR16LAW	
17" / 43cm	ALBR17LA2	ALBR17LA4	ALBR17LAR	ALBR17LAW	
18" / 45.7cm	ALBR18LA2	ALBR18LA4	ALBR18LAR	ALBR18LAW	
19" / 48.2cm	ALBR19LA2	ALBR19LA4	ALBR19LAR	ALBR19LAW	
20" / 51cm	ALBR20LA2	ALBR20LA4	ALBR20LAR	ALBR20LAW	

#### Basic Mounting

Basic Release	MSRP
ALBR10LAB	\$610.00
ALBR11LAB	
ALBR12LAB	
ALBR13LAB	
ALBR14LAB	
ALBR15LAB	
ALBR16LAB	
ALBR17LAB	
ALBR18LAB	
ALBR19LAB	
ALBR20LAB	

### Deep Contour Depth(4"/100mm), Coded E2615

#### Performance Mounting Hardware Options

Width of Wheelchair	2-point PRO	4-Point PRO	Quick Release	Width Adjustable ELITE	MSRP
10" / 25.4cm	ALBR10LD2	ALBR10LD4	ALBR10LDR	ALBR10LDW	\$685.00
11" / 28cm	ALBR11LD2	ALBR11LD4	ALBR11LDR	ALBR11LDW	
12" / 30.5cm	ALBR12LD2	ALBR12LD4	ALBR12LDR	ALBR12LDW	
13" / 33cm	ALBR13LD2	ALBR13LD4	ALBR13LDR	ALBR13LDW	
14" / 35.5cm	ALBR14LD2	ALBR14LD4	ALBR14LDR	ALBR14LDW	
15" / 38cm	ALBR15LD2	ALBR15LD4	ALBR15LDR	ALBR15LDW	
16" / 40.6cm	ALBR16LD2	ALBR16LD4	ALBR16LDR	ALBR16LDW	
17" / 43cm	ALBR17LD2	ALBR17LD4	ALBR17LDR	ALBR17LDW	
18" / 45.7cm	ALBR18LD2	ALBR18LD4	ALBR18LR	ALBR18LDW	
19" / 48.2cm	ALBR19LD2	ALBR19LD4	ALBR19LDR	ALBR19LDW	
20" / 51cm	ALBR20LD2	ALBR20LD4	ALBR20LDR	ALBR20LDW	

#### Basic Mounting

Basic Release	MSRP
ALBR10LDB	\$685.00
ALBR11LDB	
ALBR12LDB	
ALBR13LDB	
ALBR14LDB	
ALBR15LDB	
ALBR16LDB	
ALBR17LDB	
ALBR18LDB	
ALBR19LDB	
ALBR20LDB	

NOTE:

1)Each back support includes: support shell, mounting hardware, contoured foam cushion, cover, and operation manual.





User Inspired Designs

# AL Series

## Medium Height(13"/33cm), Positioning Back Supports

### Active Contour Depth(2.5"/63.5mm), Coded E2613

#### Performance Mounting

Width of Wheelchair	2-point PRO	4-Point PRO	Quick Release	Width Adjustable ELITE	MSRP
12" / 30.5cm	ALBR12MA2	ALBR12MA4	ALBR12MAR	ALBR12MAW	\$625.00
13" / 33cm	ALBR13MA2	ALBR13MA4	ALBR13MAR	ALBR13MAW	
14" / 35.5cm	ALBR14MA2	ALBR14MA4	ALBR14MAR	ALBR14MAW	
15" / 38cm	ALBR15MA2	ALBR15MA4	ALBR15MAR	ALBR15MAW	
16" / 40.6cm	ALBR16MA2	ALBR16MA4	ALBR16MAR	ALBR16MAW	
17" / 43cm	ALBR17MA2	ALBR17MA4	ALBR17MAR	ALBR17MAW	
18" / 45.7cm	ALBR18MA2	ALBR18MA4	ALBR18MAR	ALBR18MAW	
19" / 48.2cm	ALBR19MA2	ALBR19MA4	ALBR19MAR	ALBR19MAW	
20" / 51cm	ALBR20MA2	ALBR20MA4	ALBR20MAR	ALBR20MAW	

#### Basic Mounting

Basic Release	MSRP
ALBR12MAB	\$625.00
ALBR13MAB	
ALBR14MAB	
ALBR15MAB	
ALBR16MAB	
ALBR17MAB	
ALBR18MAB	
ALBR19MAB	
ALBR20MAB	

### Deep Contour Depth(4"/100mm), Coded E2615

#### Performance Mounting Hardware Options

Width of Wheelchair	2-point PRO	4-Point PRO	Quick Release	Width Adjustable ELITE	MSRP
12" / 30.5cm	ALBR12MD2	ALBR12MD4	ALBR12MDR	ALBR12MDW	\$700.00
13" / 33cm	ALBR13MD2	ALBR13MD4	ALBR13MDR	ALBR13MDW	
14" / 35.5cm	ALBR14MD2	ALBR14MD4	ALBR14MDR	ALBR14MDW	
15" / 38cm	ALBR15MD2	ALBR15MD4	ALBR15MDR	ALBR15MDW	
16" / 40.6cm	ALBR16MD2	ALBR16MD4	ALBR16MDR	ALBR16MDW	
17" / 43cm	ALBR17MD2	ALBR17MD4	ALBR17MDR	ALBR17MDW	
18" / 45.7cm	ALBR18MD2	ALBR18MD4	ALBR18MDR	ALBR18MDW	
19" / 48.2cm	ALBR19MD2	ALBR19MD4	ALBR19MDR	ALBR19MDW	
20" / 51cm	ALBR20MD2	ALBR20MD4	ALBR20MDR	ALBR20MDW	

#### Basic Mounting

Basic Release	MSRP
ALBR12MDB	\$700.00
ALBR13MDB	
ALBR14MDB	
ALBR15MDB	
ALBR16MDB	
ALBR17MDB	
ALBR18MDB	
ALBR19MDB	
ALBR20MDB	



**NOTE:**

1) Each back support includes: support shell, mounting hardware, contoured foam cushion, cover, and operation manual.



User Inspired Designs

# AL Series

## Standard Height(16"/41cm), Positioning Back Supports

### Active Contour Depth(2.5"/63.5mm), Coded E2613

#### Performance Mounting Hardware Options

#### Basic Mounting

Width of Wheelchair	2-point PRO	4-Point PRO	Quick Release	Width Adjustable ELITE	MSRP	Basic Release	MSRP
14" / 35.5cm	ALBR14SA2	ALBR14SA4	ALBR14SAR	ALBR14SAW	\$640.00	ALBR14SAB	\$640.00
15" / 38cm	ALBR15SA2	ALBR15SA4	ALBR15SAR	ALBR15SAW		ALBR15SAB	
16" / 40.6cm	ALBR16SA2	ALBR16SA4	ALBR16SAR	ALBR16SAW		ALBR16SAB	
17" / 43cm	ALBR17SA2	ALBR17SA4	ALBR17SAR	ALBR17SAW		ALBR17SAB	
18" / 45.7cm	ALBR18SA2	ALBR18SA4	ALBR18SAR	ALBR18SAW		ALBR18SAB	
19" / 48.2cm	ALBR19SA2	ALBR19SA4	ALBR19SAR	ALBR19SAW		ALBR19SAB	
20" / 51cm	ALBR20SA2	ALBR20SA4	ALBR20SAR	ALBR20SAW		ALBR20SAB	

### Deep Contour Depth(4"/100mm), Coded E2615

#### Performance Mounting Hardware Options

#### Basic Mounting

Width of Wheelchair	2-point PRO	4-Point PRO	Quick Release	Width Adjustable ELITE	MSRP	Basic Release	MSRP
14" / 35.5cm	ALBR14SD2	ALBR14SD4	ALBR14SDR	ALBR14SDW	\$715.00	ALBR14SDB	\$715.00
15" / 38cm	ALBR15SD2	ALBR15SD4	ALBR15SDR	ALBR15SDW		ALBR15SDB	
16" / 40.6cm	ALBR16SD2	ALBR16SD4	ALBR16SDR	ALBR16SDW		ALBR16SDB	
17" / 43cm	ALBR17SD2	ALBR17SD4	ALBR17SDR	ALBR17SDW		ALBR17SDB	
18" / 45.7cm	ALBR18SD2	ALBR18SD4	ALBR18SDR	ALBR18SDW		ALBR18SDB	
19" / 48.2cm	ALBR19SD2	ALBR19SD4	ALBR19SDR	ALBR19SDW		ALBR19SDB	
20" / 51cm	ALBR20SD2	ALBR20SD4	ALBR20SDR	ALBR20SDW		ALBR20SDB	

### Extra Deep Contour Depth(6"/150mm), Coded E2615

#### Performance Mounting Hardware Options

#### Basic Mounting

Width of Wheelchair	2-point PRO	4-Point PRO	Quick Release	Width Adjustable ELITE	MSRP	Basic Release	MSRP
16" / 40.6cm	ALBR16562	ALBR16564	ALBR1656R	ALBR1656W	\$790.00	ALBR1656B	\$790.00
17" / 43cm	ALBR17562	ALBR17564	ALBR1756R	ALBR1756W		ALBR1756B	
18" / 45.7cm	ALBR18562	ALBR18564	ALBR1856R	ALBR1856W		ALBR1856B	
19" / 48.2cm	ALBR19562	ALBR19564	ALBR1956R	ALBR1956W		ALBR1956B	
20" / 51cm	ALBR20562	ALBR20564	ALBR2056R	ALBR2056W		ALBR2056B	



**NOTE:**

1) Each back support includes: support shell, mounting hardware, contoured foam cushion, cover, and operation manual.



User Inspired Designs

# AL Series

## Tall Height(18"/46cm), Positioning Back Supports

### Active Contour Depth(2.5"/63.5mm), Coded E2613



Width of Wheelchair	Performance Mounting Hardware Options					Basic Mounting	
	2-point PRO	4-Point PRO	Quick Release	Width Adjustable ELITE	MSRP	Basic Release	MSRP
16" / 40.6cm	ALBR16TA2	ALBR16TA4	ALBR16TAR	ALBR16TAW	\$655.00	ALBR16TAB	\$655.00
17" / 43cm	ALBR17TA2	ALBR17TA4	ALBR17TAR	ALBR17TAW		ALBR17TAB	
18" / 45.7cm	ALBR18TA2	ALBR18TA4	ALBR18TAR	ALBR18TAW		ALBR18TAB	
19" / 48.2cm	ALBR19TA2	ALBR19TA4	ALBR19TAR	ALBR19TAW		ALBR19TAB	
20" / 51cm	ALBR20TA2	ALBR20TA4	ALBR20TAR	ALBR20TAW		ALBR20TAB	

### Deep Contour Depth (4"/100mm), Coded E2615

Width of Wheelchair	Performance Mounting Hardware Options					Basic Mounting	
	2-point PRO	4-Point PRO	Quick Release	Width Adjustable ELITE	MSRP	Basic Release	MSRP
16" / 40.6cm	ALBR16TD2	ALBR16TD4	ALBR16TDR	ALBR16TDW	\$730.00	ALBR16TDB	\$730.00
17" / 43cm	ALBR17TD2	ALBR17TD4	ALBR17TDR	ALBR17TDW		ALBR17TDB	
18" / 45.7cm	ALBR18TD2	ALBR18TD4	ALBR18TDR	ALBR18TDW		ALBR18TDB	
19" / 48.2cm	ALBR19TD2	ALBR19TD4	ALBR19TDR	ALBR19TDW		ALBR19TDB	
20" / 51cm	ALBR20TD2	ALBR20TD4	ALBR20TDR	ALBR20TDW		ALBR20TDB	

### Extra Deep Contour Depth (6"/150mm), Coded E2615

Width of Wheelchair	Performance Mounting Hardware Options					Basic Mounting	
	2-point PRO	4-Point PRO	Quick Release	Width Adjustable ELITE	MSRP	Basic Release	MSRP
16" / 40.6cm	ALBR16T62	ALBR16T64	ALBR16T6R	ALBR16T6W	\$730.00	ALBR16T6B	\$730.00
17" / 43cm	ALBR17T62	ALBR17T64	ALBR17T6R	ALBR17T6W		ALBR17T6B	
18" / 45.7cm	ALBR18T62	ALBR18T64	ALBR18T6R	ALBR18T6W		ALBR18T6B	
19" / 48.2cm	ALBR19T62	ALBR19T64	ALBR19T6R	ALBR19T6W		ALBR19T6B	
20" / 51cm	ALBR20T62	ALBR20T64	ALBR20T6R	ALBR20T6W		ALBR20T6B	



**NOTE:** 1) Each back support includes: support shell, mounting hardware, contoured foam cushion, cover, and operation manual.



User Inspired Designs

# AL Series

## Extra Tall Height(20"/51cm), positioning Back Supports

### Active Contour(2.5"/63.5mm), Coded E2613


Width of Wheelchair	Performance Mounting Hardware Options					Basic Mounting	
	2-point PRO	4-Point PRO	Quick Release	Width Adjustable ELITE	MSRP	Basic Release	MSRP
16" / 40.6cm	ALBR16XA2	ALBR16XA4	ALBR16XAR	ALBR16XAW	\$670.00	ALBR16XAB	\$670.00
17" / 43cm	ALBR17XA2	ALBR17XA4	ALBR17XAR	ALBR17XAW		ALBR17XAB	
18" / 45.7cm	ALBR18XA2	ALBR18XA4	ALBR18XAR	ALBR18XAW		ALBR18XAB	
19" / 48.2cm	ALBR19XA2	ALBR19XA4	ALBR19XAR	ALBR19XAW		ALBR19XAB	
20" / 51cm	ALBR20XA2	ALBR20XA4	ALBR20XAR	ALBR20XAW		ALBR20XAB	

### Deep Contour Depth(4"/100mm), Coded E2615

Width of Wheelchair	Performance Mounting Hardware Options					Basic Mounting	
	2-point PRO	4-Point PRO	Quick Release	Width Adjustable ELITE	MSRP	Basic Release	MSRP
16" / 40.6cm	ALBR16XD2	ALBR16XD4	ALBR16XDR	ALBR16XDW	\$745.00	ALBR16XDB	\$745.00
17" / 43cm	ALBR17XD2	ALBR17XD4	ALBR17XDR	ALBR17XDW		ALBR17XDB	
18" / 45.7cm	ALBR18XD2	ALBR18XD4	ALBR18XDR	ALBR18XDW		ALBR18XDB	
19" / 48.2cm	ALBR19XD2	ALBR19XD4	ALBR19XDR	ALBR19XDW		ALBR19XDB	
20" / 51cm	ALBR20XD2	ALBR20XD4	ALBR20XDR	ALBR20XDW		ALBR20XDB	

### Extra Deep Contour Depth(6"/150mm), Coded E2615

Width of Wheelchair	Performance Mounting Hardware Options					Basic Mounting	
	2-point PRO	4-Point PRO	Quick Release	Width Adjustable ELITE	MSRP	Basic Release	MSRP
16" / 40.6cm	ALBR16X62	ALBR16X64	ALBR16X6R	ALBR16X6W	\$820.00	ALBR16	\$820.00
17" / 43cm	ALBR17X62	ALBR17X64	ALBR17X6R	ALBR17X6W		ALBR17	
18" / 45.7cm	ALBR18X62	ALBR18X64	ALBR18X6R	ALBR18X6W		ALBR18	
19" / 48.2cm	ALBR19X62	ALBR19X64	ALBR19X6R	ALBR19X6W		ALBR19	
20" / 51cm	ALBR20X62	ALBR20X64	ALBR20X6R	ALBR20X6W		ALBR20	



**NOTE:** 1) Each Back support includes: support shell, mounting hardware, contoured foam cushion, cover, and operation manual.

# AFT Series Back Supports

**The AFT Series Back Support provides a lightweight option for those with complex positioning needs. Ideal for power or manual chair users, especially those with power tilt or tilt-in-space features.**

**Fits Chair Widths**

10" (25.4cm) - 20" (51cm)

**Support Heights**

Low: 10" (25cm)

Medium: 13" (33cm)

Standard: 16" (41cm)

Tall: 18" (43cm)

Extra Tall: 20" (51cm)

**Contour Depths**

Active: 2.5" (63.5mm)

Deep: 4" (100mm)

Extra Deep: 6" (150mm)

**Weight**

Weight: 1.5lbs (680 grams)

(approx. weight of a 16" / 40.5cm standard height back)

HCPCS Coding

Active - E2613



# The AFT Series Back Support Breakdown

Flat top design allows for greater stability throughout the scapular area when using power tilt or tilt-in-space chairs as well as help to prevent harnesses and other positioning straps from sliding down the back support

Extra tracks for additional mounting hardware adjustments, lateral support brackets, and positioning strap mounting



Slots on the top and bottom of the shell allow easy mounting of various positioning straps and 4 point harnesses

# AFT Series

## Low Height (10"/25cm), positioning Back Supports

**Deep Contour Depth (4"/100mm), Code E2615\*\***

### Performance Mounting Hardware Options

### Basic Mounting

Width of Wheelchair	Performance Mounting Hardware Options				MSRP	Basic Mounting	
	2-point PRO	4-Point PRO	Quick Release	Width Adjustable ELITE		Basic Release	MSRP
10" / 25.4cm	AFTBR10LD2	AFTBR10LD4	AFTBR10LDR	AFTBR10LDW	\$685.00	AFTBR10LDB	\$685.00
11" / 28cm	AFTBR11LD2	AFTBR11LD4	AFTBR11LDR	AFTBR11LDW		AFTBR11LDB	
12" / 30.5cm	AFTBR12LD2	AFTBR12LD4	AFTBR12LDR	AFTBR12LDW		AFTBR12LDB	
13" / 33cm	AFTBR13LD2	AFTBR13LD4	AFTBR13LDR	AFTBR13LDW		AFTBR13LDB	
14" / 35.5cm	AFTBR14LD2	AFTBR14LD4	AFTBR14LDR	AFTBR14LDW		AFTBR14LDB	
15" / 38cm	AFTBR15LD2	AFTBR15LD4	AFTBR15LDR	AFTBR15LDW		AFTBR15LDB	
16" / 40.6cm	AFTBR16LD2	AFTBR16LD4	AFTBR16LDR	AFTBR16LDW		AFTBR16LDB	
17" / 43cm	AFTBR17LD2	AFTBR17LD4	AFTBR17LDR	AFTBR17LDW		AFTBR17LDB	
18" / 45.7cm	AFTBR18LD2	AFTBR18LD4	AFTBR18LDR	AFTBR18LDW		AFTBR18LDB	
19" / 48.2cm	AFTBR19LD2	AFTBR19LD4	AFTBR19LDR	AFTBR19LDW		AFTBR19LDB	
20" / 51cm	AFTBR20LD2	AFTBR20LD4	AFTBR20LDR	AFTBR20LDW	AFTBR20LDB		



**NOTE:**

1) Each back support includes: support shell, mounting hardware, contoured foam cushion, cover, and operation manual.

2) \*\*Pending code verification.



# AL Series

## Medium Height(13"/33cm), Positioning Back Supports

**Deep Contour Depth(4"/100mm), Code E2615\*\***

### Performance Mounting Hardware Options

Width of Wheelchair	2-point PRO	4-Point PRO	Quick Release	Width Adjustable ELITE	MSRP
12" / 30.5cm	AFTBR12MD2	AFTBR12MD4	AFTBR12MDR	AFTBR12LDW	\$700.00
13" / 33cm	AFTBR13MD2	AFTBR13MD4	AFTBR13MDR	AFTBR13LDW	
14" / 35.5cm	AFTBR14MD2	AFTBR14MD4	AFTBR14MDR	AFTBR14LDW	
15" / 38cm	AFTBR15MD2	AFTBR15MD4	AFTBR15MDR	AFTBR15LDW	
16" / 40.6cm	AFTBR16MD2	AFTBR16MD4	AFTBR16MDR	AFTBR16LDW	
17" / 43cm	AFTBR17MD2	AFTBR17MD4	AFTBR17MDR	AFTBR17LDW	
18" / 45.7cm	AFTBR18MD2	AFTBR18MD4	AFTBR18MDR	AFTBR18LDW	
19" / 48.2cm	AFTBR19MD2	AFTBR19MD4	AFTBR19MDR	AFTBR19LDW	
20" / 51cm	AFTBR20MD2	AFTBR20MD4	AFTBR20MDR	AFTBR20LDW	

### Basic Mounting

Basic Release	MSRP
AFTBR12MDB	\$700.00
AFTBR13MDB	
AFTBR14MDB	
AFTBR15MDB	
AFTBR16MDB	
AFTBR17MDB	
AFTBR18MDB	
AFTBR19MDB	
AFTBR20MDB	



NOTE:

1) Each back support includes: support shell, mounting hardware, contoured foam cushion, cover, and operation manual.

2) \*\*Pending code verification.



# AFT Series

## Standard Height(16"/41cm), Positioning Back Support

### Active Contour Depth(2.5"/63.5mm), Code E2613\*\*

#### Performance Mounting Hardware Options

#### Basic Mounting

Width of Wheelchair	2-point PRO	4-Point PRO	Quick Release	Width Adjustable ELITE	MSRP	Basic Release	MSRP
14" / 35.5cm	AFTBR14SA2	AFTBR14SA4	AFTBR14SAR	AFTBR14SAW	\$640.00	AFTBR14SAB	\$640.00
15" / 38cm	AFTBR15SA2	AFTBR15SA4	AFTBR15SAR	AFTBR15SAW		AFTBR15SAB	
16" / 40.6cm	AFTBR16SA2	AFTBR16SA4	AFTBR16SAR	AFTBR16SAW		AFTBR16SAB	
17" / 43cm	AFTBR17SA2	AFTBR17SA4	AFTBR17SAR	AFTBR17SAW		AFTBR17SAB	
18" / 45.7cm	AFTBR18SA2	AFTBR18SA4	AFTBR18SAR	AFTBR18SAW		AFTBR18SAB	
19" / 48.2cm	AFTBR19SA2	AFTBR19SA4	AFTBR19SAR	AFTBR19SAW		AFTBR19SAB	
20" / 51cm	AFTBR20SA2	AFTBR20SA4	AFTBR20SAR	AFTBR20SAW	ATFBR20SAB		

### Deep Contour Depth(4"/100mm), Code E2615\*\*

#### Performance Mounting Hardware Options

#### Basic Mounting


Width of Wheelchair	2-point PRO	4-Point PRO	Quick Release	Width Adjustable ELITE	MSRP	Basic Release	MSRP
14" / 35.5cm	AFTBR14SD2	AFTBR14SD4	AFTBR14SDR	AFTBR14SDW	\$715.00	AFTBR14SDB	\$715.00
15" / 38cm	AFTBR15SD2	AFTBR15SD4	AFTBR15SDR	AFTBR15SDW		AFTBR15SDB	
16" / 40.6cm	AFTBR16SD2	AFTBR16SD4	AFTBR16SDR	AFTBR16SDW		AFTBR16SDB	
17" / 43cm	AFTBR17SD2	AFTBR17SD4	AFTBR17SDR	AFTBR17SDW		AFTBR17SDB	
18" / 45.7cm	AFTBR18SD2	AFTBR18SD4	AFTBR18SDR	AFTBR18SDW		AFTBR18SDB	
19" / 48.2cm	AFTBR19SD2	AFTBR19SD4	AFTBR19SDR	AFTBR19SDW		AFTBR19SDB	
20" / 51cm	AFTBR20SD2	AFTBR20SD4	AFTBR20SDR	AFTBR20SDW	AFTBR20SDB		

### Extra Deep Contour Depth(6"/150mm), Code E2615\*\*

#### Performance Mounting Hardware Options

#### Basic Mounting

Width of Wheelchair	2-point PRO	4-Point PRO	Quick Release	Width Adjustable ELITE	MSRP	Basic Release	MSRP
16" / 40.6cm	AFTBR16S62	AFTBR16S64	AFTBR16S6R	AFTBR16S6W	\$790.00	AFTBR16S6B	\$790.00
17" / 43cm	AFTBR17S62	AFTBR17S64	AFTBR17S6R	AFTBR17S6W		AFTBR17S6B	
18" / 45.7cm	AFTBR18S62	AFTBR18S64	AFTBR18S6R	AFTBR18S6W		AFTBR18S6B	
19" / 48.2cm	AFTBR19S62	AFTBR19S64	AFTBR19S6R	AFTBR19S6W		AFTBR19S6B	
20" / 51cm	AFTBR20S62	AFTBR20S64	AFTBR20S6R	AFTBR20S6W		AFTBR20S6B	



**NOTE:**

1) Each back support includes: support shell, mounting hardware, contoured foam cushion, cover, and operation manual.



# AFT Series

## Extra Tall Height(20"/51cm), Positioning Back Support

### Active Contour Depth(2.5"/63.5mm), Coded E2613 \*\*

#### Performance Mounting Hardware Options

#### Basic Mounting

Width of Wheelchair	2-point PRO	4-Point PRO	Quick Release	Width Adjustable ELITE	MSRP	Basic Release	MSRP
16" / 40.6cm	AFTBR16XA4	AFTBR16XA4	AFTBR16XAR	AFTBR16XAW	\$670.00	AFTBR16XAB	\$670.00
17" / 43cm	AFTBR17XA2	AFTBR17XA4	AFTBR17XAR	AFTBR17XAW		AFTBR17XAB	
18" / 45.7cm	AFTBR18XA2	AFTBR18XA4	AFTBR18XAR	AFTBR18XAW		AFTBR18XAB	
19" / 48.2cm	AFTBR19XA2	AFTBR19XA4	AFTBR19XAR	AFTBR19XAW		AFTBR19XAB	
20" / 51cm	AFTBR20XA2	AFTBR20XA4	AFTBR20XAR	AFTBR20XAW		AFTBR20XAB	

### Deep Contour Depth(4"/100mm), Code E2615\*\*

#### Performance Mounting Hardware Options

#### Basic Mounting


Width of Wheelchair	2-point PRO	4-Point PRO	Quick Release	Width Adjustable ELITE	MSRP	Basic Release	MSRP
16" / 40.6cm	AFTBR16XD2	AFTBR16XD4	AFTBR16XDR	AFTBR16XDW	\$745.00	AFTBR16XDB	\$745.00
17" / 43cm	AFTBR17XD2	AFTBR17XD4	AFTBR17XDR	AFTBR17XDW		AFTBR17XDB	
18" / 45.7cm	AFTBR18XD2	AFTBR18XD4	AFTBR18XDR	AFTBR18XDW		AFTBR18XDB	
19" / 48.2cm	AFTBR19XD2	AFTBR19XD4	AFTBR19XDR	AFTBR19XDW		AFTBR19XDB	
20" / 51cm	AFTBR20XD2	AFTBR20XD4	AFTBR20XDR	AFTBR20XDW		AFTBR20XDB	

### Extra Deep Contour Depth(6"/150mm), Code E2615\*\*

#### Performance Mounting Hardware Options

#### Basic Mounting

Width of Wheelchair	2-point PRO	4-Point PRO	Quick Release	Width Adjustable ELITE	MSRP	Basic Release	MSRP
16" / 40.6cm	AFTBR16X62	AFTBR16X64	AFTBR16X6R	AFTBR16X6W	\$820.00	AFTBR16X6B	\$820.00
17" / 43cm	AFTBR17X62	AFTBR17X64	AFTBR17X6R	AFTBR17X6W		AFTBR17X6B	
18" / 45.7cm	AFTBR18X62	AFTBR18X64	AFTBR18X6R	AFTBR18X6W		AFTBR18X6B	
19" / 48.2cm	AFTBR19X62	AFTBR19X64	AFTBR19X6R	AFTBR19X6W		AFTBR19X6B	
20" / 51cm	AFTBR20X62	AFTBR20X64	AFTBR20X6R	AFTBR20X6W		AFTBR20X6B	



NOTE:

- 1) Each back support includes: support shell, mounting hardware, contoured foam cushion, cover, and operation manual.
- 2) \*\*Pending code verification.



# **Seating**

**PSIS Support Inserts**

**Lateral Supports**

**Push Handles**

**Privacy Flaps**

**Seat Bases**

**and more...**

# Lateral Supports

Available on Active Contour AL Series and CF Series Backs\*

Integrates with mounting hardware for less weight

Custom offsets available upon request

## Mounting Options

Stationary

Swing away

\*Lateral Supports are not compatible with Fixed ELITE hardware



## Adjustable Lateral Bracket



## Pad Sizes



Small 3" H x 3" D



Medium 3" H x 4" D



Large 4.5" H x 4.5" D





X Large 4" H x 6" D

# PSIS Lumbar Support Insert

Posterior, Superior, Iliac Spine (PSIS) When positioning an ADI back, the interior portion of the back support should be positioned to make contact at the PSIS to provide optimal pelvic support and trunk control. If additional support is required, ADI offers the support insert shown below as an optional feature.

**To install, unzip the fabric cover and insert the lumbar support between the yellow foam surface and the fabric cover. Position vertically as desired.**

**Note:**  
For a secure hold, double sided Velcro is provided.





# Seating Components

## Lateral Supports for ADI Backs

(Integrates into existing mounting hardware, adjustable, active contour AL & CF back supports only)

Item #	Description	MSRP
4 LATADJSWSM	Adjustable Back - LAT Swingaway Small Pad (pad size 3" H x 3" D) (each)	<b>\$260.00</b>
4 LATADJSWMD	Adjustable Back - LAT Swingaway Medium Pad (pad size 3" H x 4" D) (each)	
4 LATADJSWLG	Adjustable Back - LAT Swingaway Large Pad (pad size 4.5" H x 4.5" D) (each)	
4 LATADJSWXLG	Adjustable Back - LAT Swingaway X-Large Pad (pad size 4" H x 6" D) (each)	
4 LATADJSTSM	Adjustable Back - LAT Stationary Small Pad (pad size 3" H x 3" D) (each)	<b>\$160.00</b>
4 LATADJSTMD	Adjustable Back - LAT Stationary Medium Pad (pad size 3" H x 4" D) (each)	
4 LATADJSTLG	Adjustable Back - LAT Stationary Large (pad size 4.5" H x 4.5" D) (each)	
4 LATADJSTXLG	Adjustable Back - LAT Stationary X-Large (pad size 4" H x 6" D) (each)	
<b>Back Support Accessories</b>		
<b>Item #</b>	<b>Description</b>	
BR6000	PSIS Support Insert, 9.5" x 3.5" x 1", each	<b>\$8.00</b>
BR6001	Loop Handle Assembly, complete	<b>\$22.50</b>
BR6002	Privacy Flap	<b>\$45.00</b>
<b>Custom Machined Push Handles</b>		
<b>Item #</b>	<b>Description</b>	
PH750	Billet push handles for 3/4" (.750") back post frame (pair)	<b>\$175.00</b>
PH875	Billet push handles for 7/8" (.8750") back post frame (pair)	
PH1	Billet push handles for 1" (1.00") back post frame (pair)	
<b>Custom Machined Key Chain Holders</b>		
KH1	Billet key chain holder 1" (1.000") frame tubing (each)	<b>\$70.00</b>
KH1125	Billet key chain holder 1 1/8" (1.125") frame tubing (each)	
KH125	Billet key chain holder 1 1/4" (1.250") frame tubing (each)	

**Carbon Fiber Seat Base and Side Guards**

**Note: CF seat base & side guards cut per chair specifications. Seat base max width 18", max depth 20"**

<b>Item #</b>	<b>Description</b>	
1	CFSB Carbon Fiber Seat Base w/ Aluminum Mounting Hardware	\$625.00
	<b>Note: mounting hardware fits round tube diameter 1" (25mm) &amp; 1.25" (32mm)</b>	
	CFSBI Carbon Fiber Seat Base Insert (no mounting hardware)	\$375.00
5	CSB Composite Seat Base, mounts to OEM seat rail insert nuts w/ under brace support	\$340.00
<b>Item #</b>	<b>Description</b>	
2	CFSGF CF Side Guard w/ folding mechanism (pair)	\$400.00
2	CFSGR CF Side Guard w/ non removable mounting hardware (pair)	\$375.00
3	CFSGS CF Side Guard to replace stock or bolt to frame / no hardware (pair)	\$250.00

**NOTE**

- 1 Each seat base includes: CF base, aluminum mounting hardware (frame tubing diameter determines mounting hardware type)
- 2 Requires CF seat base
- 3 Uses existing side guard mounting hardware
- 4 Laterals supports are universal, mounts to right or left side
- 5 Composite seat base includes: composite base, under seat support



# Disc Brakes






# Wheelchair Disc Brakes

**Please Note:**For reimbursement, use **CMS HCPCS Code: E2228** -manual wheelchair accessory, wheel braking system, complete, each

<b>DB Series: Wheelchair Disc Brake System w/ Quick Release Disc Brake Hub</b>		MSRP
<small>(System Includes: disc, caliper, caliper bracket, quick release disc hub, cable, cable housing and lever assembly)</small>		
<b>Item #:</b>	<b>Description</b>	
DBV	Disc Brake system w/ Variable Lever assembly	\$1,350.00
DBP	Disc Brake system w/ left & right side Para Levers	
\$1,150.00		
DBE	Disc Brake system w/ "E-Brakes" electric brake/lock	\$1,900.00
<b>Fusion16 DB Series Wheel</b>		
wheel features: Black AL 0M Rim, 16 Black 12G SS Spokes, AL Hub w/ integrated spline drive		
NOTE: wheel price includes - standard aluminum handrim, high-performance black/gray tire and tube		
<b>Item #:</b>	<b>Description</b>	
WF1624DB	24" - Fusion16 DB Series Wheel (each)	
WF1625DB	25" - Fusion16 DB Series Wheel (each)	
WF1626DB	26" - Fusion16 DB Series Wheel (each)	\$205.00
<b>Options</b>		
1) Hub Color: Black or Silver		
<b>Spline Drive Insert Adapter</b>		
<small>(presses into inner bearing race of wheels other than Fusion16)</small>		
<b>Item #:</b>	<b>Description</b>	
WSDI	Wheel Spline Insert - hub adapter (pair)	\$150.00



# Wheelchair Disc Brakes

<b>DM Series: Wheelchair Disc Brake System w/ Direct Mount Disc Brake Hub</b>		MSRP
(System Includes: rotor, caliper, caliper bracket, cable, cable housing, lock-up lever, wheel and solid tire)		
Item #:	Description	
DBDM12D	Disc brake / lock system for PED / Attendant chair with dual lock-up levers w/ 12" - ADI/Sun CR20 DM Series Wheels	\$750.00
DBDM16D	Disc brake / lock system for PED / Attendant chair with dual lock-up levers w/ 16" - ADI/Sun CR20 DM Series Wheels	\$765.00
DBDM20D	Disc brake / lock system for PED / Attendant chair with dual lock-up levers w/ 20" - ADI/Sun CR20 DM Series Wheels	\$780.00
DBDM24D	Disc brake / lock system for PED / Attendant chair with dual lock-up levers w/ 24" - ADI/Sun CR20 DM Series Wheels	\$795.00





User Inspired Designs

# **Wheel Components**

**Fusion16**

**Spinergy**

**Axles**

**Tires**

**Handrims**

**And more...**



# Complete Wheel Assemblies

Black Aluminum 0M Rim, 16 Black 12G Stainless Spokes, Billet Aluminum Hub		MSRP
<b>Wheel assemblies include:</b> Fusion16 wheel Standard aluminum handrim High-performance black/gray tire High-performance tube Available with black or silver hub		
<b>Item #:</b>	<b>Description</b>	
	<b>Everyday Wheel</b>	
WF1624	24" - Fusion16 Series Wheel (each)	\$205.00
WF1625	25" - Fusion16 Series Wheel (each)	
WF1626	26" - Fusion16 Series Wheel (each)	
	<b>DB Series Wheel</b>	
WF1624DB	24" - Fusion16 DB Series Wheel (each)	\$205.00
WF1625D	25" - Fusion16 DB Series Wheel (each)	
WF1626DB	26" - Fusion16 DB Series Wheel (each)	
<b>Wheel assemblies include:</b> Spinergy SPOX wheel Standard aluminum handrim High-performance black/gray tire High-performance tube		
<b>Item #:</b>	<b>Description</b>	
	<b>Spinergy 18 Spoke Wheel</b>	
WSPOX24	24" - Spinergy SPOX 18 spoke (silver hub, black spoke & rim) (each)	\$347.00
WSPOX25	25" - Spinergy SPOX 18 spoke (silver hub, black spoke & rim) (each)	
WSPOX26	26" - Spinergy SPOX 18 spoke (silver hub, black spoke & rim) (each)	
	<b>Spinergy 12 Spoke Wheel</b>	
WSPOX24LX	24" - LX Spinergy SPOX 12 spoke (silver hub, black spoke & rim) (each)	\$447.50
WSPOX25LX	25" - LX Spinergy SPOX 12 spoke (silver hub, black spoke & rim) (each)	
WSPOX26LX	26" - LX Spinergy SPOX 12 spoke (silver hub, black spoke & rim) (each)	



## Wheel Parts

Item #	Description	MSRP
WF1624B	24" - Fusion16 Wheel (black hub, black spoke & rim) - each	\$135.00
WF1624S	24" - Fusion16 Wheel (silver hub, black spoke & rim) - each	\$135.00
WF1625B	25" - Fusion16 Wheel (black hub, black spoke & rim) - each	\$135.00
WF1625S	25" - Fusion16 Wheel (silver hub, black spoke & rim) - each	\$135.00
WF1626B	26" - Fusion16 Wheel (black hub, black spoke & rim) - each	\$135.00
WF1626S	26" - Fusion16 Wheel (silver hub, black spoke & rim) - each	\$135.00
<b>Fusion16 Wheel - DB Series</b>		
WF1624BDB	24" - Fusion16 DB Wheel (black hub, black spoke & rim) - each	\$135.00
WF1624SDB	24" - Fusion16 DB Wheel (silver hub, black spoke & rim) - each	\$135.00
WF1625BDB	25" - Fusion16 DB Wheel (black hub, black spoke & rim) - each	\$135.00
WF1625SDB	25" - Fusion16 DB Wheel (silver hub, black spoke & rim) - each	\$135.00
WF1626BDB	26" - Fusion16 DB Wheel (black hub, black spoke & rim) - each	\$135.00
WF1626SDB	26" - Fusion16 DB Wheel (silver hub, black spoke & rim) - each	\$135.00
<b>Spinergy 18 Spoke Wheel</b>		
WSPOX2418	24" - Spinergy SPOX 18 spoke (silver hub, black spoke & rim) - each	\$285.00
WSPOX2518	25" - Spinergy SPOX 18 spoke (silver hub, black spoke & rim) - each	\$285.00
WSPOX2618	26" - Spinergy SPOX 18 spoke (silver hub, black spoke & rim) - each	\$285.00
<b>Spinergy 12 Spoke Wheel</b>		
WSPOX2412	24" - LX Spinergy SPOX 12 spoke (silver hub, black spoke & rim) - each	\$360.00
WSPOX2512	25" - LX Spinergy SPOX 12 spoke (silver hub, black spoke & rim) - each	\$360.00
WSPOX2612	26" - LX Spinergy SPOX 12 spoke (silver hub, black spoke & rim) - each	\$360.00
<b>Quick Release Axles</b>		
AXL004	ADI quick release axle, length 103-109.5mm (USE: Spinergy w/ DB spline drive insert only)	\$14.00
AXL005	ADI quick release axle, length 96.5-103mm (USE: Fusion DB, Fusion Everyday and Spinergy)	\$14.00
<b>KENDA - Everyday Tire</b>		
T24K184B	Kenda Speciality HP 24x1 / Iron Cap / K184 / 540 / black - each	\$28.00
T25K184B	Kenda Speciality HP 25x1 / Iron Cap / K184 / 559 / black - each	\$30.00
T26K196B	Kenda Speciality HP 26x1 / Iron Cap / K196 / 590 / black - each	\$32.00
<b>SHOX - High Pressure Solid Tire</b>		
T24SHOXB	SHOX HP - 24x1 / 540 / Black - each	\$28.00
T25SHOXB	SHOX HP - 25x1 / 559 / Black - each	\$30.00
T26SHOXB	SHOX HP - 26x1 / 590 / Black - each	\$32.00
<b>KENDA Tube</b>		
T24TSV	24 x 1" - Tube / Straight Valve / 540 - each	\$9.00
T25TSV	25 x 1" - Tube / Straight Valve / 559 - each	\$9.00
T26TSV	26 x 1" - Tube / Straight Valve / 590 - each	\$9.00



## Wheel Parts

Item #	Description	MSRP
<b>SUN Metal Handrims</b>		
HR24M	24" - Handrim, Hard Anodized Aluminum, Gray - each	\$45.00
HR25M	25" - Handrim, Hard Anodized Aluminum, Gray - each	\$45.00
HR26M	26" - Handrim, Hard Anodized Aluminum, Gray - each	\$45.00
<b>Black Vinyl Coated Handrims</b>		
HR24VC	24" - Handrim, Black Vinyl Coating, Aluminum - each	\$135.00
HR25VC	25" - Handrim, Black Vinyl Coating, Aluminum - each	\$135.00
HR26VC	26" - Handrim, Black Vinyl Coating, Aluminum - each	\$135.00
<b>Black Sparkle Foam Coated Handrims</b>		
HR24FC	24" - Handrim, Black Sparkle Foam Coating, Aluminum - each	\$135.00
HR25FC	25" - Handrim, Black Sparkle Foam Coating, Aluminum - each	\$135.00
HR26FC	26" - Handrim, Black Sparkle Foam Coating, Aluminum - each	\$135.00
<b>Push BLAX Handrims</b>		
HR24PB	24" - Handrim, BLAX Coated, Aluminum - each	\$245.00
HR25PB	25" - Handrim, BLAX Coated, Aluminum - each	\$245.00
HR26PB	26" - Handrim, BLAX Coated, Aluminum - each	\$245.00
<b>Natural-Fit Handrims (Original, Larger Profile)</b>		
HR24NF	24" - Natural Fit w/ standard grip thumb piece (oval tab) - each	\$245.00
HR25NF	25" - Natural Fit w/ standard grip thumb piece (oval tab) - each	\$245.00
HR26NF	26" - Natural Fit w/ standard grip thumb piece (oval tab) - each	\$245.00
HR24NFSG	24" - Natural Fit w/ super grip thumb piece (oval tab) - each	\$285.00
HR25NFSG	25" - Natural Fit w/ super grip thumb piece (oval tab) - each	\$285.00
HR26NFSG	26" - Natural Fit w/ super grip thumb piece (oval tab) - each	\$285.00
<b>Natural-Fit LT Handrims (Reduced Weight, Smaller Profile)</b>		
HR24NFLT	24" - Natural-Fit LT w/ standard grip thumb piece (oval tab) - each	\$265.00
HR25NFLT	25" - Natural-Fit LT w/ standard grip thumb piece (oval tab) - each	\$265.00
HR26NFLT	26" - Natural-Fit LT w/ standard grip thumb piece (oval tab) - each	\$265.00
HR24NFSGLT	24" - Natural-Fit LT w/ super grip thumb piece (oval tab) - each	\$305.00
HR25NFSGLT	25" - Natural-Fit LT w/ super grip thumb piece (oval tab) - each	\$305.00
HR26NFSGLT	26" - Natural-Fit LT w/ super grip thumb piece (oval tab) - each	\$305.00



## Wheel Parts

Item #	Description	MSRP
<b>SUN Metal Handrims</b>		
HR24M	24" - Handrim, Hard Anodized Aluminum, Gray - each	\$45.00
HR25M	25" - Handrim, Hard Anodized Aluminum, Gray - each	\$45.00
HR26M	26" - Handrim, Hard Anodized Aluminum, Gray - each	\$45.00
<b>Black Vinyl Coated Handrims</b>		
HR24VC	24" - Handrim, Black Vinyl Coating, Aluminum - each	\$135.00
HR25VC	25" - Handrim, Black Vinyl Coating, Aluminum - each	\$135.00
HR26VC	26" - Handrim, Black Vinyl Coating, Aluminum - each	\$135.00
<b>Black Sparkle Foam Coated Handrims</b>		
HR24FC	24" - Handrim, Black Sparkle Foam Coating, Aluminum - each	\$135.00
HR25FC	25" - Handrim, Black Sparkle Foam Coating, Aluminum - each	\$135.00
HR26FC	26" - Handrim, Black Sparkle Foam Coating, Aluminum - each	\$135.00
<b>Push BLAX Handrims</b>		
HR24PB	24" - Handrim, BLAX Coated, Aluminum - each	\$245.00
HR25PB	25" - Handrim, BLAX Coated, Aluminum - each	\$245.00
HR26PB	26" - Handrim, BLAX Coated, Aluminum - each	\$245.00
<b>Natural-Fit Handrims (Original, Larger Profile)</b>		
HR24NF	24" - Natural Fit w/ standard grip thumb piece (oval tab) - each	\$245.00
HR25NF	25" - Natural Fit w/ standard grip thumb piece (oval tab) - each	\$245.00
HR26NF	26" - Natural Fit w/ standard grip thumb piece (oval tab) - each	\$245.00
HR24NFSG	24" - Natural Fit w/ super grip thumb piece (oval tab) - each	\$285.00
HR25NFSG	25" - Natural Fit w/ super grip thumb piece (oval tab) - each	\$285.00
HR26NFSG	26" - Natural Fit w/ super grip thumb piece (oval tab) - each	\$285.00
<b>Natural-Fit LT Handrims (Reduced Weight, Smaller Profile)</b>		
HR24NFLT	24" - Natural-Fit LT w/ standard grip thumb piece (oval tab) - each	\$265.00
HR25NFLT	25" - Natural-Fit LT w/ standard grip thumb piece (oval tab) - each	\$265.00
HR26NFLT	26" - Natural-Fit LT w/ standard grip thumb piece (oval tab) - each	\$265.00
HR24NFSGLT	24" - Natural-Fit LT w/ super grip thumb piece (oval tab) - each	\$305.00
HR25NFSGLT	25" - Natural-Fit LT w/ super grip thumb piece (oval tab) - each	\$305.00
HR26NFSGLT	26" - Natural-Fit LT w/ super grip thumb piece (oval tab) - each	\$305.00





User Inspired Designs

# Boards & Gloves



## Anti-slip Transfer Boards

Sliding boards are used to promote independence and safety during transfers.

For those with impaired balance and upper body strength, we highly recommend ADI's **Anti-Slip Transfer Boards** as a transfer-assist tool, often providing the needed assistance to enable therapist/caregiver performance of a safer and biomechanical friendly transfer.

### Anti-Slip Material

The bottom side of each transfer board has an anti-slip surface to keep the board in place with little to no movement during each transfer.

### Hand Hole

Three options of board hand hole cutouts for optimal user manage-ability. Your choice of either ver-tical, horizontal or no hand holes.

### Hardwood Materials

Durable, layered hardwood with multiple coats of clear sealer make each board smooth and aesthetically pleasing.

### Tapered Ends

Allows the user to easily slide the board underneath his or her leg/bottom for proper placement during transfer and to aid in reduction of skin sheering



21"(53.1cm) 24"(61.0cm) 29"(73.6cm)

### Length Options



Vertical

Horizontal

None

### Hand Hole

Note: maximum weight capacity is 200lbs. (91kg)

## Push & Transfer Gloves

ADI's **Push & Transfer Gloves** have an open finger, half thumb, wrap-around design for easy don and doff. Made of leather with padded palms and wrists, the gloves save wear and tear on users' hands throughout a long day of active chair use. Three inches of Velcro closure provide a secure fit as well as allows the gloves to fit most hand sizes.



## Grip Options

All Leather / Rubber Grip

**Open finger, half thumb design enables independent donning and doffing by persons with limited hand functions**

**Made of 100% leather**

**Padded palms and wrists**

**Non-slip rubber or leather surface enhances grip for transfer, push, and braking**

**3" of dual Velcro closure for a secure fit**

Note: One size fits most



# Gloves

<b>Push &amp; Transfer Gloves</b>		
<b>NOTE: Color: Black Size: One size fits most.</b>		
<b>Item #:</b>	<b>Description</b>	<b>MSRP</b>
ADG	Wheelchair Push & Transfer Glove w/ rubber grip palm (pair) back, all leather, open finger, 2" wrist protection	
ADGL	Wheelchair Push & Transfer Glove w/ leather grip palm (pair) black, all leather, open finger, 2" wrist protection	\$60.00



# Transfer Boards


Anti-Slip Transfer Boards		MSRP
Item #:	Description	
ADV21	21" Anti-Slip transfer board w/ two vertical hand holes	\$62.00
ADV24	24" Anti-Slip transfer board w/ two vertical hand holes	
ADV29	29" Anti-Slip transfer board w/ two vertical hand holes	
ADH21	21" Anti-Slip transfer board w/ two horizontal hand holes	\$62.00
ADH24	24" Anti-Slip transfer board w/ two horizontal hand holes	
ADH29	29" Anti-Slip transfer board w/ two horizontal hand holes	
NH21	21" Anti-Slip transfer board w/ no hand holes	\$62.00
NH24	24" Anti-Slip transfer board w/ no hand holes	
NH29	29" Anti-Slip transfer board w/ no hand holes	

### Three Options In Length




### Three Options in Hand Holes

- A. Vertical
- B. Horizontal
- C. No Hand Holes



### Anatomy of an ADI Transfer Board



# Mounting Hardware

ADI offers multiple options in back support mounting hardware, each designed to fit specific needs.

Whether you want the capability to remove the backrest from the chair, need the flexibility of adjustment, or just want to "set it and forget it", we've got your mounting hardware!

## Fixed Elite™



Stationary mounting to back post, and fixed-in-place positioning (no moving parts).

### Features:

- Exclusive for Carbon Fiber (CF) Series backs
- Compatible with Active and Deep contours
- Fits on backs with 14in (35.56cm) through 18in (45.72cm) widths
- Fits on backs with 10in (25.4cm), 13in (33.02cm) and 16in (40.64cm) heights

## Width Adjustable Elite™



Stationary mounting to back post, and width adjustable fixed positioning hardware.

### Features:

- Compatible with Aluminum and Carbon Fiber backs
- Fits All Contours, All Widths, All Height

## 2-Point Pro™



Fixed mounting width adjustable positioning (angle, depth and width). Two point mounting with four point stability.

### Features:

- Performance Stationary Mounting
- Adjustable Positioning
- Angle, Width and Depth
- Compatible with all ADI Backs and Contour Depths
- 10in (25.40cm) through 23in (58.42cm) Widths and All Heights

# Mounting Hardware

## 4-Point Pro™



Fixed mounting width adjustable positioning (angle, depth and width).

Basic Mounting offers removable mounting with angle, depth and width adjustment at an economical price.

### Features:

- Performance Stationary Mounting
- Adjustable Positioning
- Angle, Width and Depth
- Compatible with All ADI Backs, Contour Depths, Widths and Heights

## Quick Release



Removable mounting width adjustable positioning (angle, depth and width).

### Features:

- Performance Mounting
- Angle: 20° of adjustment without the loss of seat depth
- Depth: Range of adjustment is 1 1/4in / 3.18cm
- Width: Approximately 1in / 2.5cm

## Basic Release



Removable mounting width adjustable positioning (angle, depth and width).

Basic Mounting offers removable mounting with angle, depth and width adjustment at an economical price.

### Features:

- Angle: 30° of adjustment without the loss of seat depth
- Depth: Range of adjustment is 1 1/4in / 5.72cm
- Width: Approximately ± 1in / 2.5cm



© 1999-2017  
1.800.965.9229 | fax 1.800.806.1225  
info@stealthproducts.com  
www.stealthproducts.com  
104 John Kelly Drive, Burnet ,TX 78611